


2020 0521 HyperMax MB+I+S+H+OS 37yoM – Effects of Hyperoxia During Exercise 2 Days in a Row

Marker	2020 0521 HyperMax Rogue Echo			2020 0522 HyperMax Rogue Echo			Reference Range
	13:30PM Pre	14:14PM After	15:15PM +30 min	9:40AM Pre	10:11AM After	10:45AM +30 min	
DNA Viruses – Plasma	ND	ND	ND	ND	ND	ND	Not Detected
DNA Parasites – Plasma	ND	ND	ND	ND	ND	ND	Not Detected
DNA Fungi – Plasma	ND	ND	3 O	1 O	ND	ND	Not Detected
DNA Bacteria – Plasma	3 O	23 O	9 O	1 O	ND	ND	Not Detected
WBC	6125	9170	7960	6858	13748	7428	4000 - 11000 Cells/uL
Lymphocyte, Absolute	1624	3114	1908	1703	5096	1730	600 - 5500 Cells/uL
Lymphocyte, Percent	27	34	24	25	37	23	10 - 45 %
CD3, Absolute (T cells)	1156	1909	1379	1178	2715	1213	606 - 3187 /uL
CD3, Percent (T cells)	68.8	60.7L	72.7	68.6	52.3L	70.2	61.0 - 80.0 %
CD4, Absolute (T helper cells)	746	1222	957	784	1520	799	365 - 2087 /uL
CD4, Percent (T helper cells)	46.0	39.2	50.1	46.0	29.8L	46.2	37.0 - 52.0 %
CD8, Absolute (T suppressor cells)	320	539	356	328	957	348	154 - 1264 /uL
CD8, Percent (T suppressor cells)	19.7	17.3	18.7	19.3	18.8	20.1	15.0 - 32.0 %
CD4/CD8 Ratio	2.3	2.3	2.7	2.4	1.6	2.3	1.0 - 3.4 Ratio
CD16/56, Absolute (NK cells)	245	719H	178	231	1777H	209	26 - 497 /uL
CD16/56, Percent (NK cells)	14.1H	22.7H	9.5	13.4H	33.6H	12.1H	3.0 - 12.0 %
CD19, Absolute (B cells)	277	449	305	282	662	257	89 - 747 /uL
CD19, Percent (B cells)	16.0	14.1	16.2	16.3	12.5	14.9	9.0 - 19.0 %
WBC	6.4	9.4	7.0	6.7	11.3H	6.8	4.0 - 11.0 k/mm3
RBC	4.86	5.22	4.87	4.62	4.87	4.62	4.30 - 6.00 m/mm3
Hemoglobin	15.2	16.4	15.3	14.4	15.0	14.1	13.0 - 18.0 g/dL
Hematocrit	44.0	48.5	43.7	41.8	45.7	41.7	40.0 - 53.0 %
MCV	90.5	92.9	89.7	90.5	93.8	90.3	78.0 - 100.0 fL
MCH	31.3	31.4	31.4	31.2	30.8	30.5	27.0 - 34.0 pg
MCHC	34.5	33.8	35.0	34.4	32.8	33.8	31.0 - 37.0 g/dL
Platelet Count	235	320	230	232	349	247	130 - 450 k/mm3
RDW(sd)	40.8	42.2	40.7	39.8	42.4	40.3	38.0 - 49.0 fL
RDW(cv)	12.4	12.4	12.4	12.2	12.2	12.2	11.0 - 15.0 %
MPV	11.7	11.7	11.4	11.3	11.7	11.6	7.5 - 14.0 fL
Segmented Neutrophils	58.7	52.9	61.3	61.9	52.6	64.2	%


Lymphocytes	27.6	34.7	25.2	25.0	35.2	23.0	%
Monocytes	9.6	9.3	9.6	9.2	9.4	9.6	%
Eosinophils	3.6	2.5	3.2	3.3	2.2	2.7	%
Basophils	0.3	0.5	0.4	0.3	0.4	0.4	%
Absolute Neutrophil	3.75	4.97	4.28	4.16	5.92	4.35	1.60 - 9.30 k/uL
Absolute Lymphocyte	1.76	3.27	1.76	1.68	3.96	1.56	0.60 - 5.50 k/uL
Absolute Monocyte	0.61	0.88	0.67	0.62	1.06	0.65	0.10 - 1.60 k/uL
Absolute Eosinophil	0.23	0.24	0.22	0.22	0.25	0.18	0.00 - 0.70 k/uL
Absolute Basophil	0.02	0.05	0.03	0.02	0.04	0.03	0.00 - 0.20 k/uL
Immature Granulocytes	0.2	0.1	0.3	0.3	0.2	0.1	%
Abs Imm Granulocytes	0.01	0.01	0.02	0.02	0.02	0.01	0.00 - 0.10 k/uL
NRBC RE, Nucleated RBC Percent	0.0	0.0	0.0	0.0	0.0	0.0	0.0 - 1.0 %
Glucose	88	119 H	104H	98	152 H	100H	65 - 99 mg/dL
Urea Nitrogen (BUN)	23	24	24	25	25	24	8 - 25 mg/dL
Creatinine	1.09	1.22	1.05	1.07	1.35	1.25	0.60 - 1.50 mg/dL
GFR Estimated	86	75	90	88	66	73	>=60 mL/min/1.73m2
BUN/Creatinine Ratio	21.1	19.7	22.9	23.4	18.5	19.2	10.0 - 28.0
Uric Acid	4.7	5.0	5.3	4.4	5.1	5.9	3.5 - 8.0 mg/dL
Sodium	140	140	138	141	143	138	134 - 147 mmol/L
Potassium	4.3	4.7	4.1	4.5	4.6	4.4	3.6 - 5.3 mmol/L
Chloride	101	100	101	104	103	100	95 - 108 mmol/L
Carbon Dioxide (CO2)	25	24	23	25	14L	21	19 - 31 mmol/L
Anion Gap	13	17	14	12	26H	17	4 - 18
Osmolality, Calculated	288	291	286	292	299H	286	275 - 295 mOsm/kg
Protein, Total	7.6	8.3H	7.5	7.2	8.0	7.2	6.0 - 8.0 g/dL
Albumin	5.2H	5.1	5.1	4.7	5.2H	4.6	3.6 - 5.1 g/dL
Globulin	2.4	3.2	2.4	2.5	2.9	2.6	1.9 - 3.7 g/dL
Albumin/Globulin Ratio	2.2	1.6	2.1	1.8	1.8	1.8	1.0 - 2.5
Cholesterol	172	197	170	157	172	154	<=199 mg/dL
Triglyceride	38	53	41	34	68	29	<=149 mg/dL
Calcium	10.5	10.8H	10.3	9.9	10.7H	9.9	8.7 - 10.4 mg/dL
Phosphorus (Inorganic)	3.7	5.1H	4.2	3.9	5.5H	3.9	2.4 - 4.8 mg/dL
Alkaline Phosphatase	40	46	40	40	48	40	40 - 140 IU/L
GGT	10	11	10	9	10	10	5 - 80 IU/L
Alanine Aminotransferase	13	18	14	15	15	14	5 - 60 IU/L

Aspartate Aminotransferase	25	23	22	19	21	22	10 - 50 IU/L
Lactate Dehydrogenase	164	196	164	138	180	159	112 - 245 IU/L
Bilirubin, Total	0.4	0.5	0.5	0.2	0.3	0.2	0.2 - 1.3 mg/dL
Cholesterol/HDL Ratio	3.0	3.2	3.0	2.9	2.9	2.9	<=4.9
HDL Cholesterol	57	62	56	55	60	54	>=40 mg/dL
Non-HDL Cholesterol	115	135H	114	102	112	100	<=129 mg/dL
LDL Cholesterol, Calculated	108 H	124H	106H	95	98	94	<=99 mg/dL
VLDL Cholesterol	8	11	8	<7	14	<7	<=29 mg/dL
VA IgM, IgG, IgA	Neg		Neg	Neg		Neg	
hs-CRP	0.4	0.5	0.4	0.4	0.5	0.4	≤0.9 mg/L
Insulin	4.2	5.9	7.6				2.6-24.9 μU/mL
Albumin	5.0	5.4H	5.0	4.8	5.4H	5.0	3.5-5.2 g/dL
Erythrocyte Sedimentation Rate	7	7	9	10	10	8	<=15 mm/hr
Estradiol	9.8L	23.8L	23.6L	12.6L	16.8L	19.1L	25.8~60.7 pg/mL
FSH	4.8	5.1	4.7	5.2	5.5	5.0	1.5~12.4 mIU/mL
DHEA-S	265.0	279.9	281.8	287.4	316.6	323.7	88.9~427.0 μg/dL
LH	2.4	2.8	2.2	4.7	4.7	3.7	1.7~8.6 mIU/mL
SHBG	77.4H	89.1H	83.6H	76.8H	88.6H	80.2H	16.5~55.9 nmol/L
Cortisol	6.6	11.6	14.6H	7.0	13.5	19.6H	6.2-19.4 μg/dL
Testosterone, Total	735.3	784.1	637.6	707.8	912.1	686.1	200.5~1437.8 ng/dL
Free Testosterone	8.46	7.91	6.72	8.25	9.52	7.57	4.09~37.37 ng/dL
Progesterone	0.083	0.246H	0.357H	0.077	0.243H	0.444H	≤0.149 ng/mL
Parathyroid Hormone	18	31	21	23	19	27	15~65 pg/mL
Estriol	<18.0	<18.0	<18.0	<18.0	<18.0	<18.0	≤200.8 pg/mL
Estrone	29.9	<9.0L	27.7	21.2	28.2	16.8	10.2~49.9 pg/mL
Prolactin	12.02	19.33H	16.30H	9.68	25.60H	21.78H	4.04~15.20 ng/mL
Dihydrotestosterone	86.5	73.2	73.0	90.2	93.4	88.6	14.8~101.8 ng/dL
Pregnenolone	0.83	1.55	0.71	0.96	0.99	0.96	0.38~3.50 ng/mL
d-ROMs <a href="https://www.hedsrl.it/eng/oxidative-stress/what-is-d-roms-test/">https://www.hedsrl.it/eng/oxidative-stress/what-is-d-roms-test/</a> Unit of measure: U. Carr 1 U. Carr = 0.08 mg H <sub>2</sub> O <sub>2</sub> /dL	414	519	464	404	548	392	250-300 Optimal value 300-320 Border line 321-340 Low ox stress 341-400 Med ox stress


							401-500 High ox stress > 500 Very high ox stress
PAT <a href="https://www.hedsrl.it/eng/oxidative-stress/what-is-the-pat-test/">https://www.hedsrl.it/eng/oxidative-stress/what-is-the-pat-test/</a> Unit of measure: U. Cor 1 U.Cor = 1.4 µMol/L of ascorbic acid	2842	2599	2971	1713	3302	1244	<2800 Very high 2200–2800 Normal 2200–2000 BL low 2000–1800 Slightly def < 1800 Deficient
OBRI <a href="https://www.hedsrl.it/eng/obri/">https://www.hedsrl.it/eng/obri/</a> Oxidative Balance Risk Index The cardiovascular risk index	1.3	1.8	1.4	2.1	1.5	2.8	0.8-1.2 Normal 1.3-1.7 Borderline 1.8-2.2 High >2.2 Very High
OSI Redox <a href="https://www.hedsrl.it/eng/osi/">https://www.hedsrl.it/eng/osi/</a> Oxidative Stress Index Summary value of oxid stress	79	139	108	98	157	120	<40 Normal 41-65 Borderline 66-120 High >121 Very High



WBC counts rise on 2<sup>nd</sup> day of hyperoxic exercise




Lymphocyte counts rise on 2<sup>nd</sup> day of hyperoxic exercise




NK Cell counts rise on 2<sup>nd</sup> day of hyperoxic exercise

### LRS




O2 saturation at 100% most of the time on both days

# LRS



EtCO2 lower on 2<sup>nd</sup> day of hyperoxic exercise



Oxidative stress hire on 2<sup>nd</sup> day of hyperoxic exercise

## Subject background:

- 37 year old male

## Exercise Protocol

- On 2020 0521
  - First blood draw at 15:30, before exercise, called Pre
  - Performed 15 min of Rogue Echo cycling exercise while breathing through mask and connected to O2 bag
  - Second blood draw at 14:14, after exercise, called After
  - Third blood draw at 15:15, 30 minutes after exercise, called +30 min
- 2020 0514
  - First blood draw at 9:40, before exercise, called Pre
  - Performed 15 min of Rogue Echo cycling exercise while breathing through mask and connected to O2 bag
  - Second blood draw at 10:11, after exercise, called After
  - Third blood draw at 11:45, 30 minutes after exercise, called +30 min

## Results:

- If you compare first day of exercise with Hypermax O2 vs second day of exercise with Hypermax O2:
  - This individual had done zero hyperoxia exercise sessions prior to this data and notice organisms were found in his blood!
  - Hypermax O2 eliminated all organisms found in plasma by second day
  - Hypermax O2 significantly increases WBCs, lymphocytes, NK cells, and platelets even more on second day!
- No negative effects of Hypermax on liver or kidney function markers, indicating this is safe
- Notice most hormones were higher on the second day of testing, especially: DHEA-S, LH, SHBG, Cortisol, Total Testosterone, Free Testosterone, Dihydrotestosterone and Pregnenolone
- Reviewing this data with respect to:
  - Immunity Boosting – data shows increase in immune parameters
  - Disease and Virus fighting/preventative – direct evidence – shows elimination of plasma organisms detected
  - Lung Health – no evidence, cannot measure adequately using Massimo
    - See <https://www.cosmed.com/en/products/cardio-pulmonary-exercise-test/quark-rmr>
    - <https://www.cosmed.com/en/products/pulmonary-function/q-box>
  - Youthfulness qualities – subjective evidence based on subject responses
    - Could be improved if we had questionnaire
    - <https://link.springer.com/article/10.1023/A:1009524612420>
    - <https://www.tandfonline.com/doi/abs/10.1080/15298868.2015.1133452?src=recsys&journalCode=psai20>
    - <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC7397859/>
  - Weight Loss – no evidence
    - We have tools to measure



- This needs to be considered as a 3-6 month project
- Reduction of Depression and Brain Fog – subjective evidence based on subject responses
  - Could be improved if we had questionnaire
  - Go here <https://www.mdcalc.com/> type in 'depression'
- Joint Pain Relief – not applicable, he has no joint issues
  - Could be improved if we had questionnaire
  - Go here <https://www.mdcalc.com/> type in 'joint pain'
- Cardiovascular Health – subjective evidence based on subject responses
  - Could be improved if we had EEG or other physiological data
- Better overall Fitness – subjective evidence based on subject responses
  - Cannot measure adequately using Massimo
  - See <https://www.cosmed.com/en/products/stress-testing-ecg>
  - <https://www.cosmed.com/en/products/ergometers/cosmed-treadmills>
  - <https://www.cosmed.com/en/products/ergometers/cosmed-bikes>